



Introducing the All-new Z-Tec Falcon XL

For Our Taller Clientele

**NEW
PRODUCT
LAUNCH**



Adjustable With Extra Height • Ultra-Lightweight • Folding



Dimensions		Features
Front Width	632mm	<ul style="list-style-type: none"> • Lightweight • Locking Handbrakes • Easy Folding • Convenient Carry Bag • Strong Build • Easy Adjusting Height • Stands When Folded • Enhanced Colours
Side Width (Wheel to Wheel)	726mm	
Handle Height (Adjustable)	780mm-932mm	
Back Width (Handle to Handle)	507mm	
Seat to Floor Height	575mm	
Seat Depth	250mm	
Weight (Including Bag)	5.9kg	

...an exciting new and enhanced colour range.



Extremely Lightweight

Our new robust design built for strength and lightweight makes transporting the Falcon XL effortless in any environment.



The Falcon XL also comes fully assembled complete with wheels attached.



Locking Handbrakes

The new Z-Tec Falcon XL features our safe locking handbrakes, providing stability for the user.

The user can squeeze the brakes to spot on the spot or push the brake handles down to lock in place.

Easy Folding

Our New Falcon XL easy folds together making it perfect for transportation.

It also stands alone when folded, so no reaching down to pick it up! And locks folded to prevent unexpected opening.



Convenient Carry Bag

Our New Falcon XL features a removable satchel bag on the front of the rollator. Capable of storing small amounts of shopping and other items.



Adjustable Height

Our New Falcon XL leads in the field of rollators for its increased height limit and the ease of adjustment.

With the height adjustment push buttons, raising and lowering a rollator has never been easier.

This is our tallest rollator yet as we have listened to our trader's concerns and the need for taller products. With simple height adjustment push buttons, raising and lowering a rollator has never been easier.



Exciting Brand New Paint Formulas

Our New Falcon XL is available in our 4 brand new colours. Rejuvenating the original colours of the original Falcon.



Fire Red



Champagne



Silver



Ice Blue

Watch this space for new colours coming soon, including our new special edition range.