

# NITRO SL ROLLATOR

---

**INNOVATIVE, STYLISH,  
COMFORTABLE AND  
SUPER LIGHTWEIGHT  
AT ONLY 5.5KG!\***



---

The super light Nitro SL was built to impress thanks to its modern design, a high-quality lightweight construction, and a range of excellent additional features as standard.

**drive**

**DeVilbiss**  
HEALTHCARE

[www.drivedevilbiss.co.uk](http://www.drivedevilbiss.co.uk)



### SUPER LIGHT

The Nitro SL is super lightweight at only 5.5kg.\*



### COMFORTABLE

The Nitro SL features Smartsoft tyres that provide a smooth and comfortable walk on a range of surfaces.



### FOLDS EASILY

The Nitro SL folds easily to a compact size for storage and includes a removable bag.



### Codes:

Code	Description
700300100	Nitro SL Rollator – Grey


### Specifications:

Length <b>72cm (28.4")</b>	Width <b>64cm (25.2")</b>	Height <b>83-93cm (32.7-36.6")</b>
Length Folded <b>72cm (28.4")</b>	Width Folded <b>27cm (10.6")</b>	Height Folded <b>83cm (32.7")</b>
Seat Height <b>57cm (22")</b>	Seat Width <b>46cm (18.1")</b>	Seat Depth <b>20cm (7.9")</b>
Width Between Handles <b>50cm (19.5")</b>	Product Weight* <b>5.5kg (12.1lb)</b>	Maximum User Weight <b>136kg (21.4st)</b>

\*Not including bag, backrest and walking stick holder. (Total weight 6.1kg)

### Features:

- » Ergonomic height adjustable handles to suit user requirements
- » Lightweight aluminium frame which has been designed for easy manoeuvrability
- » Suitable for use in tight spaces with a caster fork design that enhances turning radius
- » User friendly loop brakes with dual-function brake (for braking and parking)
- » Cross brace design for added stability
- » Easily folds with one hand to a compact size for storage and portability
- » Integral seat handle serves as a useful handgrip for pulling rollator when used as a shopper or for carrying when folded
- » Large 20cm (8") front castors for easy manoeuvrability
- » Suitable for indoor and outdoor use
- » Removable storage bag for storing shopping items and essentials
- » Plastic clip holds the rollator closed when frame is folded for easy handling
- » Backrest height is easily adjustable with tool free thumb screw
- » Soft, flexible backrest provides comfort and security while seated
- » Rear stepper tube facilitates safe curb-climbing for both attendant and user
- » Stylish tri-spoke wheel design
- » Built-in walking stick holder as standard



**If you would like to view this product please contact your local stockist.**

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.