

Dash Express

Extremely light weight, efficient opening, easy to store



Twin Brake Levers

Dual brake levers makes it easier to apply and release the brake.



Lightweight

Ultra lightweight at only 12.2kg making it easy to carry and transport.



Compact

The compact design ensures the Dash Express fits in most car boots.



DIMENSIONS	
Open width	600mm
Open length	1010mm
Open height	995mm
Seat width	17"
Folded width	260mm
Folded Length	645mm
Folded Height	765mm
Seat depth	16"
Backrest height	460mm
Backrest angle	3°
Armrest height	260mm
Front wheel	200mm
Rear wheel	300mm
Options	Dash Rehab travel bag


Weight info

Chair weight 12.2kg

Max user weight 100kg

KEY FEATURES

- ▶ Ultra lightweight
- ▶ Efficient opening
- ▶ Easy to store
- ▶ Foldable frame
- ▶ Twin brake levers



For Enquiries: