

# I Plus Power



## Specifications

FRAME	Aluminum
CONTROLLER	VR2
MOTOR	150Wx2, brushless motor
MAX SPEED	4mph (6kmh)
TRAVEL MILEAGE	9 miles (15km)
BATTERY	Lithium battery, 12Ah
FRONT WHEEL	8" x 2" PU solid castors
REAR WHEEL	12" solid PU tyres
ARMREST	Height adjustable armrest, PU armpad
FOOT REST	Detachable with angle adjustable footplate
BACK REST	Tension nylon back rest
SEAT	Breathable cushion seat
FEATURES	Safety belt; Back rest half folding

ITEM NO.	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)
SMART	100	65	95	41	45	41	48	50	120	8"	12"	23kg / 19kg

### Karma Mobility Ltd

Unit 6 Target Park, Redditch, Worcestershire B98 8YN

T: 0845 630 3436 E: info@karmamobility.co.uk

www.karmamobility.co.uk

