

## Rollators

# Days Upright Walker

Item# 7028116

The Upright Walker is a new mobility product that enables users to walk upright in a safer, more secure posture as compared to other traditional methods that cause users to hunch over. The walker is all black, made from aluminium, and supports a maximum user weight of 330 lbs (149.69 kg). Minimum assembly is required to make it as user-friendly as possible, and it comes with a cane/walking stick holder and a showerproof bag for transporting belongings.

The moulded arm rests provide extra comfort for the user, and the arms can be easily adjusted to the right length and height for each user.

The walker has smaller back wheels to improve manoeuvrability and loop lockable brakes for safety whilst walking and to ensure the walker stays in place when the user is seated. The back of the walker features two bars to assist the user in standing up.

### Key features

- Loop lockable brakes
- Comfortable seat and padded backrest
- Cane/walking stick holder
- Handles to assist standing up
- Padded forearm rests for comfort
- Adjustable height to suit the user
- Large detachable storage bag for carrying belongings
- Made from sturdy aluminium alloy
- Max user weight of 330 lbs (149.69 kg)
- Wheels are made from TPR rubber

### Dimensions

- Wheels:** 12 inch diameter, and 8 inch diameter  
(30.48 cm diameter, and 20.32 cm diameter)
- Seat size:** 18 inches wide x 11 inches deep  
(45.72 cm wide x 27.94 mm deep)
- Overall:** 33.5 inch long wheel base x 23.5 inch outside width  
(85.09 cm long wheel base x 59.69 cm outside width)



Molded arm rests for extra comfort



Cane/walking stick holder



Grip handle to aid getting in and out of the seat

