

drive

AirGo

2in1



Giving You the
Confidence to Go Further

drivedevilbiss.co.uk

AirGo 2-in-1: The Ultimate Mobility Solution

Do you love going out and about with friends and family but sometimes worry about how far you'll get on foot? With the AirGo 2-in-1, you can move confidently knowing you have a **versatile solution** at your side. This innovative product is **both a rollator and a transport chair**, allowing you to maintain independence and enjoy your day without worry. In a **matter of seconds, transition seamlessly** between **walking and sitting** in just a few **simple steps**, and rest assured that you're supported no matter where life takes you.

Open Dimensions:

Height: 78.5 - 92.5cm

Width: 67cm

Length Rollator Mode: 73cm

Length Chair Mode: 100cm

Seat Height: 51.5cm

Seat Width: 30cm

Seat Length: 50cm

Folded Dimensions:

Height: 81cm

Width: 27cm

Length: 73cm

User Weight Capacity: 136kg

Product Weight: 9.5kg*

Product Code:

700-935INT



Travel-Ready
Lightweight, foldable, and easy
to assemble without tools



5-year warranty on frame.
12 months on all other parts.
Terms and conditions apply.

drivedevilbiss.co.uk

drive



2 Step Transition
From Rollator to Transport
Chair in Seconds



Ergo-form Hand Grips
Provide easy access to brakes, with
tool-free adjustments for convenience



Patented Swing Lock Footrests
Ensuring comfort and ease of use

8" Non-marking Wheels
Perfect for indoor and outdoor use



CONTACT US FOR MORE INFORMATION
